



FAMILY
inada X.1
INADA X.1 HCP-710



Deep kneading of your neck

Firm kneading loosens your neck using INADA's unique kneading techniques.

Unique to INADA

Newly Developed


Neck Air Massage Unit

Combination technique of Neck Air Compression and Mechanical Kneading Rollers

Air bag compression loosens stiff areas around neck and shoulders.




Firm grasp



Kneading and loosening...


Squeezing the neck and applying roller massage on the top of shoulders.



Unique to INADA

Advanced stretching exercise

Whole Body Stretching



Combination of neck squeezing and whole body stretching using reclining action.

Chest and lower back Stretching



Combination of forward movement for torso and lower back with roller pressing into the center of the back.

Twisting the hips




Twisting the hips/lower body by alternating the inflation of air cells on the right and left of the seat.

Unique to INADA

Wide variety of massage techniques

Trapezius Kneading



Avoiding shoulder blades, applying deep massage into inner back muscles.

Neck Compression Massage



By removing the pillow, you can enjoy Shiatsu kneading on your neck.

Massages the very bottom of your spine



Kneading rollers reach far lower down your back with greater width of massage area. (Travel distance of 69cm with 19cm of width).

Hip / Pelvic massage



Side air cells repeat compression and release to relieve tired muscles.

Feet / Calf Massage




The calves are massaged from 3 directions (right, left and back).



Massaging the feet provides Shiatsu stimulation to the insteps and soles.

Automatic Massage Programs are available at the touch of a button.



HEALTHCARE Programs

- Full Body
- Shoulder
- Lower Back
- Stretch (Hi)
- Stretch (Lo)

QUICK Programs

- Full Body
- Shoulder
- Lower Back

INADA's unique features



Optical Shiatsu Point Sensor

Automatic Body sensing technology
Automatically detects your Shiatsu points along your back.



HEALTHCARE Programs

Professional massage techniques are pre-programmed.



Seesaw action, 4 kneading rollers & super slow movement

Soft and penetrating massage deep into the muscles.

Unique to INADA


Worlds' first development of super relaxing vibration

INADA's special vibration devices are mounted on the massage mechanism. This special feature was the result of INADA's pursuit of extreme relaxation by application of 1/f swing wave, which is unique to INADA's X.1 massage chair.

$1/f$


1/f Swing Wave

You can enjoy the gentle swing wave as if you are feeling a natural breeze or hearing the gentle murmur of a stream.




Trapezoid wave

Gentle fading out followed by gradual Vibe intensity will further enhance your relaxation.



Special kneading device with kneading rollers

The kneading rollers can deeply penetrate into your shoulders, upper back and lower back with INADA's special vibration effect.



Waveform image of Vibe 1/f



Universal Design X.1 has been designed for everybody's perfect comfort.

■ Adjustable Pillow Unit




Pillow position is fully adjustable up and down.

■ Sliding Footrest



Sliding range is 9cm, allowing anyone to fit comfortably in the leg rest.

■ Easy access



The low profile seat height has been designed for easy access.



■ Side pockets
For your convenience... not only for the remote control but magazines can be stored in the side pockets.



CREAM (CW)



BLACK (BK)

RED (RD)

Specification

Name	INADA X.1 HCP-710	
Size	Compact Size Approx. 70cm (W) x 103cm (H) x 116cm (D)	
Weight		
Power Supply	AC 220-240V	
Power Consumption	120W	
Timer	Approx. 15 min	
Reclining Angles	115°~165°	
Outer Covering	Synthetic Leather (with anti-bacterial and dirt-resistance finish)	
Accessories	Power lead, buffer pad (thick / thin), Shiatsu foot mat x 2 pcs, Shiatsu pads A and B x 2 pcs each	



Manufactured by FAMILY CO., LTD
2-1-3, Nishi-Miyahara, Yodogawa-ku, Osaka 532-0004, Japan

The duration of each massage program varies depending on height, body shape and other factors. All massage program durations are based on a person of 170 cm (5'7") height. The actual colour of the chair may differ from the printed manual. Information in this catalogue is current as of September 2009. To improve the product, specifications may be changed without notice.

Oakdale

Official European Importer/Distributor

Ellesmere Business Park
Oswestry Road, Ellesmere SY12 0EW, UK
Tel: +44 (0) 1691 623234 Fax: +44 (0) 1691 624224

Web: www.inada.co.uk Email: sales@inada.co.uk